

# REUSABLE DIAMOND BI-ASPHERIC LENSES

## OPTICS

- Patent-pending optical design
- High performance multi-layer AR and hard coatings

## ERGONOMICS

- Design patent pending
- Ultralight metallic alloy

## PERFORMANCE

- Small pupil imaging
- Low image distortion



### SMT 039-R 90D LENS

Image magnification – 0.64x  
Field of view – 75°



### SMT 042-R 28D LENS

Image magnification – 2.1x  
Field of view – 53°



### SMT 040-R 78D LENS

Image magnification – 0.77x  
Field of view – 80°



### SMT 043-R 20D LENS

Image magnification – 3.0x  
Field of view – 53°



### SMT 041-R 60D LENS

Image magnification – 0.96x  
Field of view – 68°